

Axis Scientific Classic Human Skeleton Assembly Guide

A-104267

Components:



Step 1:




Insert the black pole into the top of the base with the portion containing the black screw on top.

Step 2:




Once all 5 legs are on, flip over the base. Insert casters into the legs of the base. The casters slide easily into the hole located on the bottom of the leg.

Step 3:




Slide legs onto base piece. Be sure to slide the leg all the way in.

Step 4:



Unscrew the black screw at the top of the pole to allow the base of the skeleton to be inserted. Note: Do not attach the legs or head until the end.

Step 5:



Insert pelvic rod into the top of the stand.

Step 6:



Screw on the black screw that was previously removed. This will secure the skeleton.

Step 7:



Now attach the legs. A screw will go through the metal piece on leg so that the screw is sticking outward. Then insert the end of the screw through the hole in the pelvis.

Step 8:



Screw the nut onto the other side to secure the leg. Repeat on other side.

Step 9:



Lastly, attach the skull. There is a pin that goes through the top of the skull and can be seen inside. This must be inserted into the rod at the top of the spine.

Complete:



If all pieces are placed correctly, the skeleton will be completed and secure.

312.702.1420

AnatomyWarehouse.com

Support@AnatomyWarehouse.com